



Capsule Filling Machine NJP-400



Main Features:

1. It adopts enlarged cam, tianhong has become the one who matched the cam of RU100 to RU83, while it enhanced the stability of the machines and extends the life of the machine to a great extent.
2. It has been improved the inner design turret, and imports beeline bearing for each machine from Japan directly, so that it can be sure a longer usage life and precision of the machine.
3. It is controlled by the computer, has as steeples timing with the frequency conversion the numeral showing makes an easy operation and clear appearance.
4. The working station cam runs under the good lubricating condition, and maintains the inner lubrication of cam slot to a full extent, increasing the pressure pulverization oil pump, thus extends the operation life of the spare parts.
5. Add the extra function of capsule control switch. Control the capsule outside of the machine; it is more safely and quickly.
6. The multi-bore dosing brings a precision dosage; the good capsule applicability makes a high capsules qualified rate more than 99%. It can be filled with the Chinese traditional medicine and the western medicine.
7. Add the click-moving control system, make and easier operation and clearness.

Machinery Characteristic:

1. Working steady, safe and reliable enough.
2. Capsule qualified rate is more than 99% with the precise dosage.
3. Timing frequency-conversion, programmable controlled.
4. New closed structure, accord with GMP standard.

Technical Parameters:

1. Output: 400 capsule/min
2. No. of segment bores: 3
3. Suitable for capsule: 00-4
4. Total Powder: 3kw
5. Overall Weight: 700kg
6. Overall Dimension: 700*800*900