

GS 10-1N

Syringes filling and closing machine



Function and features:

- 1. To load syringe tubes automatically and continuously.
- 2. Watered-based solution and viscous liquid can be filled with filling space, no bubbles.
- 3. In a vacuum stoppering closed for the end of filling syringe.
- 4. Rubber plug in close contact with the outlet of the syringe without air bubbles.
- 5. Turntable structure, simple operation, automatic continuous production.

Item code	GS 10-1N
Name	Syringes filling and closing machine
Syringe volume (ml)	0.5-100
Filling volume (ml)	0-100
Filling accuracy	≤±0.3%
Speed (Pcs/minute)	30-45
Speed (Pcs/hour)	1800-2700
Total power(Kw)	2.1
Compressed air	0.55-0.75 Mpa 20L/S
Overall dimension (L*W*H)	1500*850*1800
Weight of basic machine (kg)	600