

HPT-II Capsule Polisher with Sorter



The Model HPT-II polisher with the function of sorter as well. It plays the function of not only the polishing for capsule and tablet but eliminating static electricity. It can also automatically reject low weight capsule; loose piece and fragments of capsules. It is suitable for all sizes of capsules without the need to change parts. Sorter: This capsule sorter is used to select out the empty capsules, loose pieces as well as scraps after capsule filling. It can also automatically select out those capsules with lighter dosage. The machine can be connected to the main production line. The capsules can be polished after sorting process which is more suitable for SR capsule dosages.

Feature:

- The height and angle can be adjusted to fit into any chute of the capsule filling machines and tablet presses so the products can be polished immediately after production.
- It can eliminate static; and automatically reject ultra-lightweight capsules, fragment, and empty capsules.
- The new type of net cylinder ensures no jammed capsules during operations.
- The capsules are not directly contacted with the metal net to protect the printed capsule effectively.
- The new type of brush is durable and can be changed easily.
- Excellent design for quick cleaning and straightforward maintenance.
- Adopts frequency converter, which is great for continuous long hours of operations.
- Drive by synchronous belt to reduce noise and vibration of the machine.
- It is suitable for all sizes of capsules without any change parts.
- All the main parts are made of premium stainless steel are in compliance with pharmaceutical GMP requirements.

Technical Parameters:

Output	Max 7000 (Doc/min)
Power	0.2kw
Power	220v 50Hz
Overall dimensions	1300x500x1200 mm
Weight	50Kg