

## *GMP Oral Liquid Filling Line*



This GMP oral liquid filling line is specially designed for washing, high-temperature sterilization drying, filling, and sealing 10ml-30ml oral liquid bottles. This line is in accordance with GMP standards, and is quite suitable for measured filling for oral liquid bottles.

Its main parts are ultrasonic washer, high temperature sterilization tunnel oven, and filling and sealing machine. And its different machines can be used individually, and they should adopt box packaging and container transport.

### **Features:**

1. Controlled by computer and touch screen, this production line enjoys high automation, and its output can be as many as 36000 bottles per hour.
2. Easy to install, and convenient to maintain.

### **Main Parameters:**

|                     |                           |
|---------------------|---------------------------|
| Bottle Size         | 10ml~30ml straight bottle |
| Production Capacity | 6000-36000 bottle/hour    |
| Power               | ≤60KW                     |
| Electric Supply     | 380V/50hz, 3phase         |
| Size                | 9600×1400×2000mm          |