

## *ZPW-21K Rotary Tablet Press*



This machine is used to press all sorts of granular materials, like geometry shape of profiled sheet, double lettering piece, traditional Chinese medicine tablet, western Medicine tablet, salt tablet, calcium tablet, buccal tablet, spirulina tablet, milk tablet, candy tablet, effervescent tablets, chewable tablets, seasoning block, chicken pieces, camphor ball(naphthalene), aromatic piece, disinfection tablet, fertilizer tablet, pesticide tablet, button cell, electronics, metallurgy powder ceramic grain etc.

### **Features**

- The structure has double pressure; two sets of roller, frequency control of motor speed, depth of fill, tablet thickness can be adjusted.
- It equipped with mechanical buffer device, the pressure overload can automatically shut down, and killing all had a broken impact in accordance with GMP.
- According to user's requirements, the machine can be equipped with a vacuum cleaner, screen machine, feeding machine etc.

## Technical Parameters

Model	ZP21K
Dies	21 sets
Max. Pressure	120 KN
Max. Diameter of Tablet	30 - 40 mm
Max .Depth of Fill	35 mm
Max .Thickness of Tablet	17 mm
Turret Speed	20 r/min
Max .Production Capacity	50400 tablets/hour
Motor	4 KW 380V/50HZ 3 phase
Overall Size	1200 x 1200 x 1800 mm
Net Weight	2600 kg