

The Dividella *NeoTOP 104*



Range of formats:

min: 60x45x17 – max.: 240x170x70

Pack contents / range of volumes:

Blisters

Parenterals:

Ampoules / vials cartridges / syringes
bottles

Medical devices: inhalers, nasalsprays, pens,
etc.

Output:

up to 10 Cartons

Packaging:

Printed or plain carton blanks from flat die cut
blanks. Partition walls from flat die cut blanks.

Format change:

10-20 min - (1 Person)
no tools required

Application:

Semi automatic forming and erecting of
NeoTOP cartons including integrated partition
from flat blanks.

Cylindrical objects, applicators, tubes, medical
devices or blisters are then laid in manually.
At the closing station the packs are tucked in
and sealed by glueing station.

The machine is designed for applications
where flexibility for small lots or clinical trials
are required.

Machine size:

Carton Erecting Module:

approx. 2300 x 1750 x 1900 mm

**Transfer belt for manual insertion of
objects:**

Depending on requirements

Carton Closing Module:

approx. 1960 x 1100 x 1120

