

Nail Polish Filling And Capping Machine XQZGX-JY12



Application:

- 1. It is widely used in precision, quantitative loading / dilution of pharmaceutical, pharmaceutical, food, chemical, cosmetic and other industries.
- 2. Pre-chromatography preparation, pharmaceutical QC and R&D, environmental monitoring, mineral and blood analysis, forensic identification.
- 3. Used for precise quantitative loading / dilution of blood, biological reagents, oral agents, injections, advanced cosmetics, flavors, fungicides and other liquids.

Characteristic

- 1. Intelligent operating system with touch screen is used.
- 2. Liquid filling part: constant pressure feeding, suitable for filling medium and high viscosity liquids.
- 3. 304 stainless steel appearance, desktop structure, practical and convenient.
- 4. The part of the capping is turned down by pneumatic torque to ensure the consistency of the tightness of the cover.
- 5. This machine adopts the way of quick filling of bottle molds, which can be easily and quickly replaced among different bottle types.
- 6. This equipment adopts German sensor to realize the functions of no filling without bottle, no capping without lids.
- 7. Can be customized according to the need to increase the plug or cover and other functions.
- 8. Taiwan brand pneumatic components are adopted.

Technical parameters

Model	XQZGX-JY12
Manufacturing speed	0- 35PCS / minute
Tank capacity	SS304 5L / 10L / 20L
Size of Host Machine	500*550*650(mm)
Power	AC 220V 50-60Hz 450W
Air pressure	0.4-0.6Mpa