

SED-100RJ-B Fish Oil Soft Gel Capsule Encapsulation Machine For Any Laboratory



The SED-100RJ-B soft capsule machine is a new type of product carefully designed by our company based on the soft capsule production process and production characteristics in accordance with GMP requirements. The whole equipment adopts optimized design and all uses 220V power supply, which is simple to operate. It is mainly used for the development of soft capsule products in laboratories and pharmaceutical research units or small batch production trial production.

Advantages

- 1) Shrinkable design according to soft capsule production process, complete functions.
- 2) Compact structure, the entire equipment covers an area of only 1.5m².
- 3) The whole set of equipment can be easily moved and is suitable for any laboratory.

Product Parameters

Model	SED-100RJ-B
Rolling mold size	64*100mm
Rolling mold speed	0~4rpm
Feed amount	6 *(0-1.5) mL
Total Weight	400kg
Dimensions	620 *660 *1350mm
Type	Soft Gelatin Encapsulation Machine