

SED-110RJ Laboratory Use Small Softgel Encapsulation Machine for Small Scale Production



Soft gel encapsulation machine can quantificational inject the oil, mixing suspension and pasty mass into the capsule to form all kinds of soft gel encapsulation with different size, shape and color.

Because of the merits of high speed disintegration, easy to oral, hard to oxygenation, easy to store and take, the soft gel encapsulation is welcomed by more and more customers.

Owing to the high productivity, charging exactly, quality stabilization, high finished product rate, it is favored by more and more manufacturers and used in the fields of medicine, health products, cosmetic and game-color ball and so on.

This product has passed provincial-level appraisal of new product and the certification of EU adopted two items of patent technology in manufacturing with exquisite workmanship.

Feature:

1. Optimal design, compact structure, and small footprint.
2. Meet the GMP standard.
3. All using a touch screen PLC control, convenient operation, and stable performance.
4. The entire device can be easily moved.
5. The gelatin box adopts liquid sensor, control the quantity of the liquid gelatin automatically.
6. The whole machine adopts newest temperature controlling system, so the temperature can be controlled.

Technical Data

Model	SED-110RJ
Rotary speed of die roll	0-7rpm
Supplying Volume for flow single-pillar pistol	0-1ml
Difference among filling volume	±2%
Power for whole machine	2.3KW
Pillar quantity	8
Main Body sizes	1680×1100×1820
Weight of whole machine	600Kg