

## *Blister Packing Machine BM-III*



This machine is designed for small pharm factory, healthcare plant, food plants and preparation rooms, can be used to pack solid drugs like capsule, tablet, gelatin pearl, suppository, milk tablet, candy and hardware. The whole set is of compact design, beautiful appearance, smooth operation and low noise. Use internal heating method to ensure the PVC is evenly heated. There is no bubble in the aluminum foil while heating and hence no waste will occur.

### **Specifications**

<b>Model</b>	<b>BM-III</b>
<b>Max. Heating Area(mm)</b>	80*100
<b>Overall Dimensions(mm)</b>	250*200*235
<b>Power(W)</b>	500
<b>Capacity(pcs/min)</b>	5--20
<b>Package Size(mm)</b>	370*330*310
<b>Weight(kg)</b>	14