



## ***VPF-TCI Automatic Vial Powder Filling & Stoppering Isolator (Weighing System)***



### **Specification**

Vial Size: Ø 22 mm ~ Ø 30 mm; height 35 mm ~ 75 mm

Stopper Size: Ø 20 mm

Capacity: 60 ~ 100 pcs/min

(Depends on filling volume, product characteristics & specific weight, vial outer diameter and height)

Filling Range: 40 mg ~ 2000 mg

(4 different ranges of dosing disc between 40 mg ~ 2000 mg for filling; equip with multiple filling system, range can be increased to higher weight.)

The filling dosage depends on product characteristics and specific weight; client shall supply powder sample for testing.

Filling Accuracy: ±3% (Based on 1 g and should depends on filling dosage, product characteristics and specific weight.)