



## ***RJWJ-PG soft capsule polishing equipment***



RJWJ-PG soft capsule polishing equipment for polishing the post-forming treatment, the treated surface gloss bright soft capsule, pill and beautiful, and easy to use and intestinal digestion and absorption.

### **Technical Parameters:**

<b>Model</b>	<b>PG1000</b>	<b>PG800</b>
The diameter of the pot body :	1000mm	800mm
Producing capacity :	50-80kg/time	30-60 kg/time
The tilt angle $\alpha$ :	20-45 <sup>0</sup>	20-45 <sup>0</sup>
Rotating speed of the spot :	26round/min	30round/min
Dimension :	1100x1100x1800mm	925x925x1600mm
Weight :	260kg	250kg