



Tea Bag Packing Machine CYBZ-01



The tea bag packing machine is designed to pack tea bags with tea in them. The machine can finish bag making, measuring, filling, sealing, threading feeding, tag cutting and counting etc. It is suitable for package of broken tea, coffee, ginseng essence, diet tea, health-caring tea, medicine tea, as well as tea leaves and herb beverage.

Features:

1. Adopt PLC to control the machine, convenient to operate.
2. The parts connected with products are made of SUS304 stainless steel, which can ensure the products healthful and reliable.
3. The machine can automatically finish bag making, measuring, filling, sealing, thread feeding, tag cutting and counting.

Introduction:

A tea bag is a small, porous sealed bag containing tea leaves and used for brewing tea, or herbs for brewing tisanes. The bag contains the tea leaves while the tea is brewed, making it easier to dispose of the leaves, and performs the same function as a tea infuser. Some tea bags have an attached piece of string with a paper label at the top that assists in removing the bag while also identifying the variety of tea. The tea bags are usually square or rectangular in shape.

Application:

Drinking tea has great health benefits. Drinking tea with tea bag is more convenient. Also tea bag has many other benefits.

1. If your skin has burned by the sun. Wet tea bags applied to the affected areas will soothe the skin.
2. Soak two tea bags in lukewarm water and place them over your closed eyes for 20 minutes. The tannins in the tea act to reduce puffiness and soothe the tired eyes.
3. Wash your hair with a cup of freshly brewed tea as a final rinse after shampooing. This will give natural shine to the hair and them silky.
4. Boil three or four tea bags in one quart of water for ten minutes. Let it cool so that it makes a comfortable soak. Soak your feet for twenty minutes. And it kills the bacteria which cause the foot odor.
5. Soak a tea bag and apply directly over the canker sore for several minutes. The tannin helps healing the sore.

Technical Parameters:

Model	Capacity	Power Source	Weight	Dimensions
CYBZ-01	28-50 bags/min	220V/1.8Kw	400kg	1050*910*2010(mm)