

SHH-GYT8 Oral Liquid Filling Capping Machine



It's just for 10-20ML oral filling&sealing, a major machine in production line for oral liquid, it main used in filling, capping placing and sealing of straight tube bottle. Tracking filling, cap automatically conveying by electromagnetic oscillator, three knives centrifugal sealing (or screwing), no bottle no filling. Filling and sealing integrated in one machine, advanced design, compact structure and completely comply with requirement of GMP.this kind of machine.

This machine is mainly used for syrup, drinks, soybean sauce, vinegar, cod-liver oil, olive oil, liquid calcium, wine, liquid fertilizer, essential oil, toner, hair oil, laundry detergent, ink, disinfectant, injection, industrial alcohol and various medical liquid filling and capping.

Characteristic:

- ◆All elements is well know brand.
- ◆Direct viewing operation, no bottle no filling, measuring accurate.
- ◆Three knives centrifugal sealing, ensure firm and beauty.
- ◆The whole machine is designed according to the requirements of GMP.

Features:

1. Adopts plunger type metering pump filling, high precision; The structure of pump adopts fast disassembly institutions, easy to clean and disinfect.
2. PLC control system, color touch screen display, machine has function of counting automatically.
3. The filling nozzle was equipped with drip-proof device.

4. When filling liquid, the filling nozzle dive into the bottom of bottle, rising slowly, it can effectively prevent bubbles.
5. Machine body is made of 304 stainless steel, full compliance with GMP standard.

Technical Parameters:

Model	SHH-HYT8
Filling volume	10-20ml
Capacity	100-150BPM
Filling Precision	≥99%
Power supply	380V, 50HZ
Sealing rate	≥99%
Power	2.0KW
Weight	750KG
Dimension	2500*1000*1800