



## **KGS10**

Antibiotic vial liquid filling and stopper-adding(full and half) linkage line

## Product Elucidation:

This linkage line consists of KCQ60 type vertical ultrasonic washing machine, GMS-B type hot air circulatory tunnel sterilizing drier and KGS10 antibiotic vial liquid filling and stopper-adding(full or half) machine. This line is widely used for freeze-drying product and blood product manufacturers, and it fully comply with GMP requirements.

Main Technical Specifications	
Applied bottle range	2-100ml
Filling error	≤2%
Whole line bottle breakage ratio	≤1%
Auto cap-adding ratio	≥ 99.5%
Whole line noise	≤70db
Powersupply	380V 50HZ
Total power	69kw
Dimension	9100*1900*2150 mm
weight:	6800kg