



GL2-25 Laboratory Dryer Granulator



Dryer granulator suitable for laboratory, middle laboratory tester and small batch manufacture, it just need 100grams powder into making forming, get request granule, you can adjust granule size and compactness, by PLC control auto control, it will suitable difference requirements, it will improve finish rate, high efficiency, high automatic, easy to use and maintenance, lower noise. Suitable pharmaceutical and food and use daily goods

This machine meet with GMP manufacture made, it is SUS304 stainless steel contact with product, pressure roller by special wear-resisting material, structure simple easy to clean, convenience to out feeder, no buildup dead

This machine meet with GMP manufacture made, it is SUS304 stainless steel contact with product, pressure roller by special wear-resisting material, structure simple easy to clean, convenience to out feeder, no buildup dead

1. This machine by PLC process control, with alarm system, avoid press roller locked rotor, alarm and removal of fault
2. Store information menu in order to control difference material craft data
3. Manual and auto operate

Main technical:

Item	unit	GL2-25
Roller size	$\Phi \times W$ mm	80×25
Capacity:	kg/h	1~5
Speed	rpm	0~18
Screw rotate speed	rpm	0~61
Granule rotate speed	rpm	0~600
Finish size	mm	0.3~1.5
Roller pressure	MPa	0~7
Total power	kw	2.63
Outside size	L×W×H mm	750×650×1350
Weight	kg	200