

## ZPW-25K Special Salt Rotary Tablet Press



This machine is used to press all sorts of granular materials, like geometry shape of profiled sheet, double lettering piece, traditional Chinese medicine tablet, western medicine tablet, salt tablet, calcium tablet, buccal tablet, spirulina tablet, milk tablet, candy tablet, effervescent tablets, chewable tablets, seasoning block, chicken pieces, camphor ball(naphthalene), aromatic piece, disinfection tablet, fertilizer tablet, pesticide tablet, button cell, electronics, metallurgy powder ceramic grain etc.

### Features

- The structure has double pressure; two sets of roller, frequency control of motor speed, depth of fill, tablet thickness can be adjusted.
- It equipped with mechanical buffer device, the pressure overload can automatically shut down, and killing all had a broken impact in accordance with GMP.
- According to user's requirements, the machine can be equipped with a vacuum cleaner, screen machine, feeding machine etc.

### Technical Parameters

| Model                    | ZPW-25K                |
|--------------------------|------------------------|
| Dies                     | 25 sets                |
| Max. Pressure            | 120 KN                 |
| Max. Diameter of Tablet  | 30 mm                  |
| Max .Depth of Fill       | 36 mm                  |
| Max .Thickness of Tablet | 15 mm                  |
| Turret Speed             | 20 r/min               |
| Max .Production Capacity | 60000 tablets/hour     |
| Motor                    | 4 KW 380V/50HZ 3 phase |
| Overall Size             | 1200 x 1200 x 1800 mm  |
| Net Weight               | 2600 kg                |